


Matrix Material

MATRIX BANDS AND RETAINERS

PATTERSON



Tofflemire Matrix Bands Stainless Steel

Precontoured bands are easier to form for excellent proximal contacts. Mirror finish for easy removal.

12/Package	1@ \$1.20
#1, .0015 Gauge088-6614	10@ \$1.05
#2, .0015 Gauge088-6622	
#3, .0015 Gauge088-6630	
#13, .0015 Gauge088-6648	
#1, .002 Gauge.....088-6655	
#2, .002 Gauge.....088-6663	
#3, .002 Gauge.....088-6671	
#13, .002 Gauge088-6689	



Tofflemire Matrix Retainers Stainless Steel

Stainless steel retainer is easy to adjust and provides a positive grip. Package contains retainer and assortment of 1 dozen bands.

Package	1@ \$14.80
Universal.....088-6515	3@ \$14.25
Contra Angle088-6523	
Junior Contra Angle..088-6531	

AutoMatrix® - Dentsply/Caulk Retainerless Matrix System

For use in posterior restorations.

Introductory Package

Contains: 24 medium regular bands, 24 narrow regular bands, 24 wide regular bands, 1 Automate Dual Torque Tightening Device and 1 snippers.

Pkg. 168-7680 **1@ \$520.00**

AutoMatrix Band Refill **1@ \$70.50**

72/Pkg. Medium Thin **3@ \$68.25**

1/4" x .0015" 168-8209

Narrow Regular

3/16" x .002" 168-8233

Medium Regular

1/4" x .002" 168-8266

Wide Regular

3/16" x .002" 168-8290

AutoMatrix® - Dentsply/Caulk Clear Retainerless Matrix Accessories

AUTOMATE® III Tightening Device

Each 168-5155 **1@ \$349.95**

Snippers

Each 168-9769 **1@ \$170.00**

Extra Shaft Assembly

Each 168-9777 **1@ \$94.00**

Contoured Metal Matrix Bands - Premier

30/Pkg
0.035mm570-9399 **1@ \$29.90**

Cure-Thru® Matrix Bands - Premier

Mylar
100/Package **1@ \$107.75**
Premolar580-6823
Molar580-6815



Composit-Tight® 3D - Garrison Dental

Sectional matrix systems are unbeatable for predictable contacts on posterior composites, and Composit-Tight sets the standard. NOW, with Composit-Tight 3D, you can achieve easy and fast tight, anatomical contacts with almost-zero learning curve! The bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour. The band thickness is .0013" (30 µm)

3D Kit with forceps (3DK2)

Kit326-6947 **1@ \$409.20**

Contains: 3 Soft Face Rings, 2 thin tine G-Rings, 50 bicuspid bands, 25 extended small bands, 50 mid molar bands, 50 molar bands, 25 extended molar bands and a specially designed Ring Placement Forceps.

3D Kit

Kit (3DK1)326-6970 **1@ \$409.20**

Contains: 3 Soft Face Rings, 2 thin tine G-Rings, 50 bicuspid bands, 25 extended small bands, 50 mid molar bands, 50 molar bands, 25 extended molar bands

Matrix bands refills

M050 3.2mm Pedodontic with extension
50/Pkg.....326-7127 **1@ \$68.35**

M100 3.8mm Bicuspid
100/Pkg.....326-7135 **1@ \$68.35**

M150 3.8mm Bicuspid with extension
50/Pkg.....326-7143 **1@ \$68.35**

M175 5.5mm Molar
100/Pkg.....326-7150 **1@ \$68.35**

M200 6.4mm Molar
50/Pkg.....326-7168 **1@ \$68.35**

M300 6.4mm Molar with extension
50/Pkg.....326-7176 **1@ \$68.35**

Soft-Face 3D-Ring (3D500)

2/Pkg.....326-6996 **1@ \$136.15**

3D Thin Tine G-Ring (3D400)

3/Pkg.....326-6988 **1@ \$85.80**

Ring Placement Forceps (AUMRDF-100)

Each326-6152 **1@ \$80.85**

Getz® Contour Matrix Bands – Waterpik

Pre-contoured band made of high tensile strength stainless steel. Will not tear under stress.

50/Pkg

.002 Regular.....336-7315	1@ \$56.85
.0015 Ultrathin336-7323	




ORIGINAL Tofflemire II Retainers – Waterpik

Polished Stainless Steel

- Tighten and release easily
- Clean, smooth surface is free of burs and sharp edges
- Contra-Angle model can be used both buccally and lingually

Package contains retainer only, no bands.
 Universal.....336-5921 1@ \$31.25
 Contra Angle336-5954



Contact Matrix Introductory Kit – Danville

Clinical Kit236-3588 1@ \$146.55

Contains: 3 each inward and 3 outward rings, 100 small, 100 large thin-flex matrices.

Masters Kit236-3570 1@ \$281.10

Contains: Clinical Kit, 1 ring pliers and 1 ring forceps.

Contact Ring Refill

Contains: 1 inward/1 outward.

Each236-3620 1@ \$27.65

Thin-Flex Matrices

100/Pkg

Small236-2218	1@ \$51.45
Large.....236-2226	3@ \$49.50



Palodent® Sectional Matrix System – Dentsply Caulk

Sectional Matrix System

For gentle separation during cavity preparation.

Palodent® Masters Kit

Contains: 100 standard Palodent® Matrices, 50 Mini-Matrices, 50 Palodent® Plus-Matrices, 4 Original BiTine® Rings & 2 BiTine II® Rings. Instructions included. Also available with a nine-minute demonstration video.
 Kit w/o video549-3721 1@ \$174.00

Palodent® Matrix Refills

Standard size Palodent® Matrices for most posterior restorations.

100/package.....549-3655 1@ \$55.15

Palodent® Mini-matrices

A smaller matrix for limited exposures.

50/package.....549-3713 1@ \$37.40

Palodent® Plus Matrix

Large matrices for subgingival or broadly extended restorations.

50/package.....549-3754 1@ \$45.95
 3@ \$45.05

BiTine® Ring Refill

4/pkg.....549-3663 1@ \$36.25

BiTine II® Ring Refill

3/pkg.....549-3697 6@ \$28.25

Tofflemire Matrix Bands – Waterpik

Original "Dead Soft" .001 gauge.

36/Package 1@ \$11.00

#2 Adult336-8289

#1 Adult336-8297

#13 Pedo336-8305

Ho Bands – Young

.001 Thin

Tofflemire Matrix bands.

Regular, 100/Box

#1472-2716 1@ \$42.80

#2472-2724

#13472-2732

"Dead Soft"

#1472-2781 1@ \$47.55

#2472-2799

#13472-2807



OptraMatrix – Ivoclar Vivadent

Selectively thinned matrix bands



Refills, Pack of 501 @ \$40.30

Molar

(Mod-Univ.) 714-1526

(I + III + D/II + IV M)....714-1534

(I + III + M/II + IV D)....714-1542

Pre-molar

(Mod-Univ.) 714-1559

(I + III + D/II + IV M)....714-1567

(I + III + M/II + IV D)....714-1575

ORIGINAL® Tofflemire Matrix Bands – Waterpik

Stainless Steel

36/Package 1@ \$6.75

#1 Universal

.0015ga.....336-6028

#2 MOD Wide

.0015ga.....336-6044

#3 MOD Narrow

.0015ga.....336-6069

#13 Pedo Universal

.0015ga.....336-6093

#14 Pedo MOD Narrow

.0015ga.....336-6127

#15 Pedo MOD Wide

.0015ga.....336-6135

#1 Universal

.002ga.....336-6010

#2 MOD Wide

.002ga.....336-6036

#13 Pedo Universal

.002ga.....336-6085



Supermat® Matrix System – Kerr

Designed for posterior macro restorations.

A ring matrix system incorporating both steel and transparent plastic. Handy tensioning device is small in size allowing for multiple simultaneous placements.

Contains: 1 SuperLock® tensioning instrument, 10 SuperCap® steel matrices, 10 SuperCap® transparent plastic matrices.

Kit433-3902 1@ \$255.45

50/Pkg

SuperCap® Matrice Refill

5mm Metal433-3928 1@ \$72.70

6.3mm Metal433-3936 2@ \$70.30

5mm Clear.....433-3944

6.3 mm Clear.....433-3951

100/Pkg

SuperCap® Spools

5.9mm Blue.....433-4512 1@ \$45.50

6.7mm Green433-4520 2@ \$44.20

MATRIX BANDS AND RETAINERS (CONT.)



"T"-Bands – Pulpdent

Self-contained matrix bands that do not require matrix retainers. Available straight, curved, narrow, wide and assorted in either brass or stainless steel.

Box/100 .002" thick 1@ \$17.10
Curved

Assorted Stainless.....591-3553
Assorted Brass.....591-3736
Narrow Stainless.....591-3439
Narrow Brass.....591-3611
Wide Stainless.....591-3496
Wide Brass.....591-3678

Straight 1@ \$17.10

Assorted Stainless.....591-3520
Assorted Brass.....591-3702
Narrow Stainless.....591-3405
Narrow Brass.....591-3587
Wide Stainless.....591-3462
Wide Brass.....591-3645

Stiff Flex Matrices – Danville



100/Pkg
Small236-2234 1@ \$51.45
Large.....236-2242 3@ \$49.50

T-Bands – P.N. Condit

Matrix Bands

Regular 100/Box 1@ \$39.60
Straight Large218-4612
Straight Medium218-4646
Straight Small218-4679
Straight Assorted.....218-4703
Curved Small218-4737
Curved Large218-4760
Curved Assorted218-4794



STRIPS AND FORMS

PATTERSON



Mylar® Matrix Strips

4" x 3/8" x .002"

Made of DuPont Mylar®.

100/Pkg.....089-7157 1@ \$4.30

6@ \$3.65

500/Pkg.....089-7165 1@ \$10.15

6@ \$9.30



Matrix Strip Holder Clip

Stainless Steel

45°.....089-7116 1@ \$18.85

90°.....089-7124 3@ \$17.30

Accor Matrix System – Premier

Multi-Sized Core Matrix Forms

Can be used on anterior and posterior teeth, with glass ionomers, composites, pattern acrylic or softened wax. Disposable.

50/Pkg.....570-9357 1@ \$46.15

Contour™ Strips – Ivoclar Vivadent

Mylar® Strip

A specially designed Mylar® Strip to fit into the sulcus, retracting the gingiva. Also helps to restrict contaminating gingival cervical fluid.

100/Pkg.....708-6093 1@ \$27.80

3@ \$26.45



Cure-Thru® Clear Cervical Matrices – Premier

Light-Cure Matrix

Compress composite material during light cure Class V procedures. Excellent for use with glass ionomers.

150/Pkg

719 Ct575-9196 1@ \$114.10

720 Ct575-9204 3@ \$112.30

721 Ct575-9212

722 Ct575-9220

723 Ct575-9238

Assorted Package

Contains: 250 matrices, 50 each of 5 sizes above.

Pkg.....575-9154 1@ \$175.60



Cure-Thru® Squeeze Matrix – Premier

Retainerless Matrix

Thin and anatomical – .075mm. They are easy to place and have a soft metal ring which, when squeezed, tightens the clear band to the tooth, which eliminates the need for matrix retainers.

Intro Package

Contains: 50 Cure-Thru Squeeze Matrices (25 each molar and premolar), and 50 Cure-Thru Wedges (25 each small and medium).

Pkg.....575-9097 1@ \$200.30

Band Refill 1@ \$132.40

50/Package 2@ \$130.30

Premolar.....575-9113

Molar575-9121

Epitex™ Matrix Strip Refill – GC America, Inc.

Clear

Contains: 1-10 meter reel.

Pkg.....207-2973 1@ \$37.10

Matrix Dispenser – Zirc

Allows matrix to be contoured as it dispenses. Contains one roll of 3/8" x 12 yd matrix tape, 0.002 gauge stainless steel cutter.
Blue.....771-6483



1@ \$23.00

Polyester Matrix Tape – 3M ESPE

With dispenser
Pkg., #1939.....501-6241

1@ \$55.00
3@ \$54.00



Stop Strips – Premier

Mylar® Anterior Bands with Integrated Stopper

100/Pkg
10mm.....576-4428 1@ \$35.00
8mm.....576-4410

Assorted Package

Contains: 75 - 10mm and 75 - 8mm.
Pkg.....576-4477 1@ \$52.50

WEDGES

Chevron Wedges – Dux

Plastic, hollow, V-shaped wedges allow for easy insertion.

400/Box
Small
Translucent Blue.....162-7595 1@ \$89.10
Large
Clear162-7603 1@ \$87.25



Contact Wedge – Danville

Provides complete gingival margin adaptation and seal. Simultaneously wedges from both buccal and lingual. Is compatible with all interproximal matrix systems.
Available in 3 sizes.

Kit
Contains: 85 of each size
Kit234-3804 1@ \$131.95

Refills of 85 1@ \$50.10
Small 2.0mm234-3820
Medium 2.6mm234-3838
Large 3.2mm234-3846

Cure-Thru Wedge System – Premier

Transparent Wedges

Intro System

Contains: 50 each of small and medium wedges, 30 each of premolar and molar bands.
Pkg.....586-1901 1@ \$169.95
3@ \$165.00



Cure-Thru Reflective Curing Wedges – Premier

Reflective Interior

White plastic interior deflects 90% of light at a 90° angle into the proximal surface for a thoroughly cured restoration without marginal gaps.

100/Pkg
Small586-1927 1@ \$104.65
Medium.....586-1935 3@ \$103.30
Ultra-Thin.....586-1943



Sycamore Wedges – Premier

Wood Wedges

Anatomically carved. Wedge number designates length of wedge in mm.

400/Pkg
#13XT (Neutral)586-0416 1@ \$50.85
#13T (Green)586-0341 3@ \$49.50
#13 (Yellow)586-0358
#15 (Blue)586-0366
#16 (Red)586-0374
#17 (Neutral).....586-0382
Assorted586-0390



Wedge Organizer – Premier

Each586-1976 1@ \$23.60
*Wedges not included.

Wizard® Anatomical Wedges – Waterpik

Wood Wedges

Non-splintering birchwood. Colour-coded for convenience.

400/Package 1@ \$52.95
3@ \$51.60
Small (Natural).....341-9611
Medium (Yellow).....341-9629
Large (Purple)341-9637
X-Large (Red)341-9645
Assorted341-9652



-  Small
-  Large
-  Wide Base
-  Slim Jim

Wizard® Wedges – Waterpik

Wood Wedges

Made of white birchwood. Triangular in cross-section, longitudinally tapered.

500/Package 1@ \$57.00
3@ \$54.80
Small342-0528
Large.....342-0551
Slim Jim342-0585
Wide Base342-0619
Assorted342-0643