

# Miscellaneous



## LIGHT BULBS

### Operatory Light Replacement Lamp Chart

| Mfg.      | Model                  | Mfg. Code       | EFOS# | Patterson# | Price      |
|-----------|------------------------|-----------------|-------|------------|------------|
| Belmont   |                        | H3              | 4071  | 70130023   | 1@ \$7.00  |
| Castle    | Panovision             | 125T10P         | 4072  | 70130032   | 1@ \$26.40 |
| Chayes    | Virginia               | H3              | 4071  | 70130023   | 1@ \$7.00  |
| Ellman    | Conversion             | Q250/MC/CL/120v | 4073  | 70130068   | 1@ \$33.60 |
| Eurodent  |                        | FCS             | 4075  | 70130024   | 1@ \$7.75  |
| F.A.R.O.  | Koolbeam               | DZE/FDS         | 4074  | 70130022   | 1@ \$21.45 |
| F.A.R.O.  | Sunlight 70 (Upgraded) | DZE/FDS         | 4074  | 70130022   | 1@ \$21.45 |
| F.A.R.O.  | Sunlight 90            | DZE/FDS         | 4074  | 70130022   | 1@ \$21.45 |
| Fiad      | The Light              | H3              | 4071  | 70130023   | 1@ \$7.00  |
| Healthco  | Aristocrat             | H3              | 4071  | 70130023   | 1@ \$7.00  |
| MDT       | Rolux                  | 12-2024-10      | 3017  | 70130028   | 1@ \$32.45 |
| J. Morita |                        | FCS             | 4075  | 70130024   | 1@ \$7.75  |
| Pelton    | Light Fantastic        | Q150T4/CL 25v   | 4076  | 70130027   | 1@ \$34.15 |
| Porter    | Daystar S-90B          | DZE/FDS         | 4074  | 70130022   | 1@ \$21.45 |
| Ritter    | J Starlite             | Q75T3/CL/DC 28v | 4079  | 70130021   | 1@ \$27.90 |
| Ritter    | K Starlite             | DZE/FDS         | 4074  | 70130022   | 1@ \$21.45 |
| Weber     | Skylight SK900G        | EPV Sk580       | 4080  | 70130052   | 1@ \$31.55 |
| Weber     | Skylight SK100         | FDT 12v 100w    | 4081  | 70130020   | 1@ \$29.80 |
| Barnstead | EJL F/Rolux Light      | 12202410        | 3017  | 70130028   | 1@ \$32.45 |
| Dental Ez | Daray                  | DRK1024         |       | 70130033   | 1@ \$93.00 |
| Pelton    | F/Executive F13 "D"    |                 | 3052  | 70130073   | 1@ \$20.90 |
| Sirona    | Bulb for C8 Unit       | 4889965         |       | 56070415   | 1@ \$16.05 |

### Light Reflectors - PAW

#### For Pelton & Crane LFI

Each .....56010086 1@ \$459.50

#### For Pelton & Crane LFIII

Each .....56010084 1@ \$502.80

#### For Ritter

Each .....56010098 1@ \$643.25

#### For Marus

Each .....56010090 1@ \$643.25

### Replacement Lamp Chart

| Mfg.        | Model                        | Mfg. Code      | EFOS# | Patterson# | Price       |
|-------------|------------------------------|----------------|-------|------------|-------------|
| 3M ESPE     | Ortholux                     | 78-8045-4571-9 | 3016  | 70130049   | 1@ \$45.90  |
| 3M ESPE     | Visilux 1/2/<br>XC 3000/2500 | 78-8045-4571-9 | 3016  | 70130049   | 1@ \$45.90  |
| 3M ESPE     | Visilux 1/2/<br>XC 3000/2500 | 78-8045-4571-9 | 3016  | 70130049   | 1@ \$45.90  |
| 3M ESPE     | XL1500 Lamp<br>52W 10V       | 9726-3         |       | 70130165   | 1@ \$40.40  |
| 3M ESPE     | F/Visio Elipar               | 77160          |       | 70130046   | 1@ \$189.45 |
| A-Dec       | Bulb for #898                | 0.022076       |       | 70130037   | 1@ \$89.00  |
| A-Dec       | F/695/6/896/8/TA-96/8L       | 0.022076       |       | 70130037   | 1@ \$89.00  |
| A-Dec       | F/Lt25 Low Speed Motor       | 0.019055       |       | 70130215   | 1@ \$103.40 |
| A-Dec       | F-F/O HPC<br>(1988-Present)  | 041.317.00     |       | 56070078   | 1@ \$107.00 |
| A-Dec       | F-F/O HPC<br>(Before 1988)   | 041.117.00     |       | 56070054   | 1@ \$107.85 |
| A-Dec       | PHIL. Quartz Halo HIX 14623  | 041.709.00     |       | 70130451   | 1@ \$34.10  |
| Air Tech    | Vista Cam - 150W             | 77121          |       | 70130232   | 1@ \$86.45  |
| Almore      | OP II                        | 3442           | 3016  | 70130049   | 1@ \$45.90  |
| Belmont     | Halogen H3<br>12V 55W        | 800-01006      |       | 70130023   | 1@ \$7.00   |
| Bosworth    | Superlite                    | 92102          | 3014  | 70130389   | 1@ \$33.90  |
| Caulk       | Max-Lite                     | 644725         | 3090  | 70130422   | 1@ \$38.95  |
| Caulk       | Prisma-Lite                  | 644101         | 3017  | 70130028   | 1@ \$32.45  |
| Caulk       | Prismetics-Lite              | 644358         | 3021  | 70130406   | 1@ \$66.65  |
| Caulk       | ProLite                      | 644815         | 3029  | 70130275   | 1@ \$42.70  |
| Caulk       | QHL75/<br>Spectrum 800 75W   | 644657         |       | 70130258   | 1@ \$125.25 |
| Caulk       | Spectrum                     | 644751         | 3009  | 70130165   | 1@ \$40.40  |
| Caulk       | Triad VLC                    | 95775          |       | 70130042   | 1@ \$50.35  |
| Coe         | Coe-Lite                     | 4070110        | 3014  | 70130389   | 1@ \$33.90  |
| Clive Craig | Viewer                       | F14T8D         |       | 70130059   | 1@ \$18.50  |
| Darby       | Superdent 2001               | 950-9023       | 3016  | 70130049   | 1@ \$45.90  |
| Darby       | Superdent 7500               | 950-9024       | 3022  | 70130174   | 1@ \$38.20  |
| Demetron    | Optilux 101                  | 20437          | 3016  | 70130049   | 1@ \$45.90  |
| Demetron    | Optilux 400                  | 20437          | 3016  | 70130049   | 1@ \$45.90  |

**Replacement Lamp Chart (Cont.)**

| Mfg.         | Model                                | Mfg. Code  | EFOS# | Patterson# | Price       |
|--------------|--------------------------------------|------------|-------|------------|-------------|
| Demetron     | Optilux 401, 501                     | 20437      | 3016  | 70130049   | 1@ \$45.90  |
| Demetron     | Optilux High Output<br>401, 403, 501 | 20140      | 3021  | 70130406   | 1@ \$66.65  |
| Demetron     | Optilux 50                           | 20919      | 3026  | 70130275   | 1@ \$42.70  |
| Demetron     | VCL 200                              | 20437      | 3016  | 70130049   | 1@ \$45.90  |
| Demetron     | VCL 300                              | 20437      | 3016  | 70130049   | 1@ \$45.90  |
| Denmat       | Ceri-Light                           | 3961       | 3021  | 70130406   | 1@ \$66.65  |
| Denmat       | Visar II                             | 3904       | 3022  | 70130174   | 1@ \$38.20  |
| Dental Comp  | Enz Lamp<br>for Rolux II             | 12-1909-20 | 4084  | 70130072   | 1@ \$52.30  |
| Dentamerica  | Litex<br>660, 680                    | 668        | 3016  | 70130049   | 1@ \$45.90  |
| Dolan-Jenner | Fiberlite 170                        | EJV        | 3022  | 70130174   | 1@ \$38.20  |
| EFOS         | EFOS LITE                            | 3021       | 3021  | 70130406   | 1@ \$66.65  |
| EFOS         | Penmed 200L<br>Orig Wand             | 640-30280  | 3022  | 70130174   | 1@ \$38.20  |
| Getz         | Focus Activator                      | 4501       | 3021  | 70130406   | 1@ \$66.65  |
| Healthco     | Aristocrat                           | HCOB400491 | 3016  | 70130049   | 1@ \$45.90  |
| Healthco     | Aristocrat II                        | HCOB400351 | 3026  | 70130275   | 1@ \$42.70  |
| Healthco     | VLC Illuminator                      | HCOB400481 | 3015  | 70130029   | 1@ \$44.05  |
| Kavo         | Light Source                         |            | 3015  | 70130029   | 1@ \$44.05  |
| Kavo         | Lux Coupler 2/Pkg                    | 553-3881   | 3090  | 70130422   | 1@ \$38.95  |
| Kavo         | Vicon DLS                            | 555-5001   | 3015  | 70130029   | 1@ \$44.05  |
| Kavo         | Vicon LS                             | 555-5001   | 3015  | 70130029   | 1@ \$44.05  |
| Kavo         | Vicon Polar II                       | 555-5001   | 3015  | 70130029   | 1@ \$44.05  |
| Kerr         | Command Light                        | 15512      | 3011  | 70130030   | 1@ \$30.80  |
| Kinetic      | Sunlite II                           | 1450       | 3022  | 70130174   | 1@ \$38.25  |
| Kinetic      | Sunspot                              | 1511       | 3021  | 70130406   | 1@ \$66.65  |
| Kulzer       | Translux CL                          | 2244       | 3021  | 70130406   | 1@ \$66.65  |
| Kulzer       | Translux II                          | 2045       | 3014  | 70130389   | 1@ \$33.90  |
| Lares        | Quartz Fire (SX)                     | 10037      | 3015  | 70130029   | 1@ \$44.05  |
| Litema       | Pluraflex HL-150                     | EFR        | 3014  | 70130389   | 1@ \$33.90  |
| Midwest      | F/5 Hole & 6 Pin                     | 790257     |       | 70130293   | 1@ \$144.10 |
| Midwest      | F/P Optic<br>w/QD Hose               | 650064     |       | 70130088   | 1@ \$123.25 |
| Midwest      | F/P Optic<br>w/WO Swivel             | 650197     |       | 70130087   | 1@ \$119.30 |
| Midwest      | Insight 1                            | 490874     | 3015  | 70130029   | 1@ \$44.05  |

| Mfg.        | Model                       | Mfg. Code   | EFOS# | Patterson# | Price       |
|-------------|-----------------------------|-------------|-------|------------|-------------|
| Midwest     | Insight 2                   | 490874      | 3015  | 70130029   | 1@ \$44.05  |
| Midwest     | Insight I                   | 490874      | 3015  | 70130029   | 1@ \$44.05  |
| Midwest     | Insight II                  | 490874      | 3015  | 70130029   | 1@ \$44.05  |
| Midwest     | Oral Illuminator 2          | 490709      | 3024  | 70130039   | 1@ \$59.45  |
| Midwest     | Oral Illuminator II         | 490709      | 3024  | 70130039   | 1@ \$59.45  |
| Midwest     | Oral Illuminator 3          | 490709      | 3024  | 70130039   | 1@ \$59.45  |
| Midwest     | Oral Illuminator III        | 490709      | 3024  | 70130039   | 1@ \$59.45  |
| Miltex NSK  | Opticure                    | EPF/64627   | 3020  | 70130170   | 1@ \$54.15  |
| O'Ryan      | Omega 110                   | ONE         | 3019  | 70130060   | 1@ \$60.65  |
| O'Ryan      | Omega 510                   | ONE         | 3019  | 70130060   | 1@ \$60.65  |
| Patterson   | TCL-390                     | 640-3026P   | 3026  | 70130275   | 1@ \$42.70  |
| Patterson   | TCL-490                     | 921311      |       | 70130257   | 1@ \$54.05  |
| Penn-Med    | Curing Gun 31L              | 301315      | 3016  | 70130049   | 1@ \$45.90  |
| Penn-Med    | Curing Wand<br>200L         | 301203      | 3022  | 70130174   | 1@ \$38.20  |
| Phasealloy  | Witelite                    | 3051        | 3028  | 70130077   | 1@ \$54.85  |
| Pro-Den     | Executor                    | 20-055      | 3021  | 70130406   | 1@ \$66.65  |
| Proma       | Power Optic Bulb            | 92475K      |       | 70130178   | 1@ \$150.30 |
| S.S.White   | Profile VCS                 | 17244       | 3021  | 70130406   | 1@ \$66.65  |
| Schein      | VCL 300                     | 100-5657    | 3021  | 70130406   | 1@ \$66.65  |
| SciCan      | Dual Lix 3 (Curing)         | EJA         | 3018  | 70130014   | 1@ \$41.35  |
| Solid State | Initiator                   | 6220002-290 | 3014  | 70130389   | 1@ \$33.90  |
| Star        | 6 Pin Swivel                | 62522       |       | 70130254   | 1@ \$155.05 |
| Star        | F/Fiber Optic "L"           | 54536       |       | 60070527   | 1@ \$168.00 |
| Star        | KRPTN/Hal Fbr Optc          | 60986       |       | 60070527   | 1@ \$168.00 |
| Star X-Ray  | Viewer                      | F1478D      |       | 70130059   | 1@ \$18.50  |
| Vivadent    | Astralix<br>Halogen 75W 12V | 551730      |       | 70130289   | 1@ \$74.20  |
| Vivadent    | F/Single<br>Function Lite   | 6010126     |       | 70130014   | 1@ \$41.35  |
| Vivadent    | Heliolux 75W 12V            | 134171      |       | 70130289   | 1@ \$74.20  |
| Vivadent    | Heliolux Lite               | 133762      | 3026  | 70130275   | 1@ \$42.70  |
| Vivadent    | Heliomat (Multi)            | 7110        | 3014  | 70130389   | 1@ \$33.90  |
| Vivadent    | Heliomat (Single)           | 7210        | 3018  | 70130014   | 1@ \$41.35  |
| Vivadent    | Vivalux VLC<br>35W 6V       | 105851      | 3037  | 70130410   | 1@ \$75.85  |
| Welch       | Denlite                     | DP1020      |       | 70130205   | 1@ \$42.10  |
| Whaledent   | Coltolux<br>3, 4 & 75       | C7925       | 3016  | 70130049   | 1@ \$45.90  |

**AIR FRESHENERS**

**Airiell™ Air Freshener  
- Ecolab**

Combines the benefits of an attractive, compact dispenser with the time-release technology of Airiell refills. Together, they provide fresh fragrances and odour counteraction for a wide variety of areas and spaces.

Dispenser size: 3¾" x 4¼"

Refills: 6/Pkg

1@ \$77.95

- Vanilla .....376-2069
- Encounter .....376-2051
- Spice Air .....376-2119
- Sea Air .....376-2127



**DAPPEN DISHES**

**PATTERSON**



**Dappen Dishes  
Glass**

- Each
- Amber.....084-9141
- Blue.....084-9158
- Green.....084-9166
- Clear .....084-9174

1@ \$3.60  
6@ \$3.35

**Dapaway  
Disposable Dishes  
- Dentsply**

**Dappen Dish  
Refill Package**

Contains: 1,000 dappen dishes, 1 finger holder.  
Pkg.....231-6198 1@ \$82.25



**Dispensing Cups -  
Young**

Plastic calibrated cups (30ml) for dispensing powders, liquids, and gels in the laboratory or at chairside. Calibrated in ml, cc, oz, and drams.  
100/Pkg.....472-0157 1@ \$16.00



# Miscellaneous

## MEDICAMENTS



### Emergency kit - BDM

|                        |          |             |
|------------------------|----------|-------------|
| Kit .....              | 133-0869 | 1@ \$219.95 |
| Refill kit .....       | 133-0851 | 1@ \$179.95 |
| Ammonia Inhalants...   | 133-1156 | 1@ \$9.60   |
| Atropine Sulfate ..... | 133-1057 | 1@ \$37.50  |
| Diphenhydramine .....  | 133-0794 | 1@ \$40.40  |
| Epinephrine .....      | 133-1032 | 1@ \$38.95  |
| Insta-Glucose .....    | 150-2533 | 1@ \$19.40  |



### Ora5® - McHenry Laboratories

Topical Bactericidal Agent For Oral Mucosa

|                 |          |            |
|-----------------|----------|------------|
| 1oz Bottle..... | 481-7953 | 1@ \$37.30 |
|                 |          | 3@ \$33.95 |

### Patient Size Bottles

|                      |          |            |
|----------------------|----------|------------|
| 1/2oz Bottles        |          |            |
| 10 Bottles/Pkg ..... | 481-7995 | 1@ \$83.25 |

### Ora5 Placement Instrument - Premier

|            |          |            |
|------------|----------|------------|
| Each ..... | 580-8357 | 1@ \$36.25 |
|------------|----------|------------|

## NITROUS

### Breathing Bag - Porter

|                      |          |            |
|----------------------|----------|------------|
| 3.5 Litre/Each ..... | 71010038 | 1@ \$63.95 |
|----------------------|----------|------------|



### DynoMite Nasal Hood - Matrix

#### Disposable

- Reduces risk of cross infection
- X-ray compatible
- Designed for a soft snug fit to put patients at ease

#### Small Size; 24/Pkg

|                  |          |            |
|------------------|----------|------------|
| Bubble Gum ..... | 480-9810 | 1@ \$99.90 |
| Strawberry ..... | 480-9828 |            |
| Orange .....     | 480-9836 |            |
| Vanilla.....     | 480-9844 |            |
| Plain .....      | 480-9851 |            |
| Assorted .....   | 480-9869 |            |

#### Medium Size; 24/Pkg

|                  |          |            |
|------------------|----------|------------|
| Bubble Gum ..... | 480-9877 | 1@ \$99.90 |
| Strawberry ..... | 480-9885 |            |
| Orange .....     | 480-9893 |            |
| Vanilla.....     | 480-9901 |            |
| Plain .....      | 480-9919 |            |
| Assorted .....   | 480-9927 |            |

#### Large Size; 24/Pkg

|                  |          |            |
|------------------|----------|------------|
| Bubble Gum ..... | 480-9935 | 1@ \$99.90 |
| Strawberry ..... | 480-9943 |            |
| Orange .....     | 480-9950 |            |
| Vanilla.....     | 480-9968 |            |
| Plain .....      | 480-9976 |            |
| Assorted .....   | 480-9984 |            |



### Hood Liner - Porter

|             |          |            |
|-------------|----------|------------|
| Adult Size  |          |            |
| 3/Pkg ..... | 71010121 | 1@ \$86.85 |

#### Pedo Size

|             |          |            |
|-------------|----------|------------|
| 3/Pkg ..... | 71010220 | 1@ \$86.85 |
|-------------|----------|------------|

### Inhaler Liners - White Gaps, Inc.

Disposable infection control barriers for nitrous oxide hoods.

|              |          |            |
|--------------|----------|------------|
| Adult size   |          |            |
| 50/Box ..... | 758-3156 | 1@ \$75.00 |
| Pedo size    |          |            |
| 50/Box ..... | 758-3164 | 1@ \$75.75 |



### Comfort Cushion - Mission Dental

#### Semi-Disposable

A uniquely designed nasal analgesia circuit. Sets new standards in comfort, convenience, cleanliness, and efficiency. A lightweight plastic, semi-disposable nasal hood, delivery and scavenging system for nitrous oxide.

|             |           |            |
|-------------|-----------|------------|
| Stac-Pac    | 3/Package | 1@ \$64.05 |
| Child ..... | 509-7621  |            |
| Adult.....  | 509-7613  |            |

|                 |          |            |
|-----------------|----------|------------|
| Individual Mask |          | 1@ \$56.25 |
| Adult without   |          | 3@ \$55.00 |
| Scavenger.....  | 509-7217 |            |

|                 |          |           |
|-----------------|----------|-----------|
| Head Strap only |          | 1@ \$9.20 |
| 2/Pkg .....     | 509-7266 |           |

|                          |          |            |
|--------------------------|----------|------------|
| Scavenger Conversion Kit |          | 1@ \$13.70 |
| Each .....               | 509-7209 |            |



### Nasal Hood - Porter

For Porter/Brown Scavenger System

|                 |          |             |
|-----------------|----------|-------------|
| Pedo Hood.....  | 71010074 | 1@ \$311.10 |
| Adult Hood..... | 71010073 | 1@ \$311.10 |



**Personal Inhaler Plus – Accutron  
Single Use Nasal Hoods**

Reassure your patients by protecting them from the risks of cross infection. Save your staff the time and effort involved in cleaning and sterilizing expensive multi-use nasal hoods. Available in a variety of flavours and sizes.

**Small Size 24/Package 1@ \$113.75**

- Groovy Grape.....102-2342
- Outlaw Orange.....102-2292
- Sassy Strawberry.....102-2300
- Birthday Bubblegum..102-2318
- Unscented
- Whisper Grey.....102-2334
- Variety Pack 1.....102-2326  
(Orange, Strawberry, Bubblegum)

**Medium Size 24/Package 1@ \$126.90**

- Outlaw Orange.....102-2110
- Sassy Strawberry.....102-2128
- Birthday Bubblegum..102-2136
- Pina Colada.....102-2144
- French Vanilla.....102-2151
- Fresh Mint.....102-2193
- Perfect Peach.....102-2201
- Unscented
- Whisper Grey.....102-2185
- Variety Pack 1.....102-2169  
(Outlaw Orange, Sassy Strawberry, Birthday Bubblegum)
- Variety Pack 2.....102-2177  
(Pina Colada, French Vanilla, Outlaw Orange)
- Variety Pack 3.....102-2219  
(Fresh Mint, Perfect Peach, French Vanilla)

**Large Size 24/Package 1@ \$133.45**

- Fresh Mint.....102-2375
- Perfect Peach.....102-2383
- French Vanilla.....102-2391
- Unscented
- Whisper Grey.....102-2409
- Variety Pack 3.....102-2417  
(Fresh Mint, Perfect Peach, French Vanilla)

**PATIENT COMFORT**



**Memory Backrest™ – Crescent Products**

Great for supporting the lumbar area of the back during lengthy procedures.

- Teal .....232-1487 1@ \$70.35
- Beige .....232-1495 1@ \$70.35



**Memory Headrest™ – Crescent Products**

Shapes to the head and neck for optimum patient comfort. Designed to fit and attach quickly to all major dental chairs.

- Teal .....232-1503 1@ \$77.30
- Beige .....232-1511 1@ \$77.30



**Memory Bodyrest™ – Crescent Products**

- Beige.....232-1537 1@ \$260.00
- Teal .....232-1529 1@ \$248.75

**TOYS**



**Itty Bittys – Rand**

*Rubber Animal Toys*

- Assorted  
50/Pkg.....592-4998 1@ \$30.00



**“My Tooth” Boxes – Coralite**  
*For Children Teeth, 1 1/8” Square*

- 100/Package 1@ \$47.50
- Clear .....224-8524
- Clear / Blue.....224-8616

**Tooth Fairy Bank – The Perfect Present**



Today's tender moments can be tomorrow's treasured memories.

- Size: 6" x 6" x 4"
- Blue.....569-0011 1@ \$27.20
- Pink.....569-0029
- Yellow.....569-0037
- Mint.....569-0045



**Treasure Chest Toys – Rand**  
*Treasure Chest*

- With 200 Toys 1@ \$85.25
- Pkg.....592-5425
- Toy Assortment Treasure Chest 1@ \$79.95
- Refill 200/Pkg.....592-5433 2@ \$78.70

## MISCELLANEOUS

### E-Z Access Shelf – Zirc

- Provides E-Z access to work space for auxiliary items
- E-Z mounting to standard equipment and adjustment to multiple positions



- Weight capacity of 25lb
  - Stain resistant finish; surface disinfectants only
- |            |          |             |
|------------|----------|-------------|
| Grey ..... | 771-0726 | 1@ \$132.00 |
| Beige..... | 771-0718 |             |
| White..... | 770-9017 |             |

### Shelf Adapter

Fits posts 1 3/8" x 1 1/4".

- |            |          |           |
|------------|----------|-----------|
| Grey ..... | 771-0742 | 1@ \$4.50 |
| Beige..... | 771-0759 |           |
| White..... | 771-0767 |           |



### Surface Mount

Mount to horizontal flat surfaces only.

- |            |          |            |
|------------|----------|------------|
| Beige..... | 771-2672 | 1@ \$17.00 |
| Grey ..... | 771-2680 |            |
| White..... | 771-2664 |            |



### 8" Extension Arm

- May be used with wall mount, surface mount or directly on to post
- Shelf has strength to hold 25lb when connected to arm
- Enables shelf to extend over patient

- |            |          |            |
|------------|----------|------------|
| Beige..... | 770-9074 | 1@ \$34.00 |
| Grey ..... | 770-9082 |            |
| White..... | 770-9066 |            |



### Wall Mount

- Allows you to mount the shelf onto a wall or side of the cabinet
- 90° angle of mount permits 180° movement of shelf on cabinet or next to wall

- |            |          |            |
|------------|----------|------------|
| Beige..... | 770-9132 | 1@ \$34.00 |
| Grey ..... | 770-9140 |            |
| White..... | 770-9124 |            |

- |                          |          |            |
|--------------------------|----------|------------|
| Clip Only.....           | 770-9165 | 1@ \$16.35 |
| HVE Holder .....         | 770-9181 | 1@ \$28.00 |
| Saliva Ejector Holder... | 770-9173 | 1@ \$28.00 |
| Cup Holder.....          | 771-0783 | 1@ \$11.50 |



### AED Plus – Zoll

*The first and only full-rescue AED*

The Zoll AED Plus™ does more than just defibrillate. It is the only AED that insures infrequent rescuers respond appropriately during a sudden cardiac arrest emergency. How? It coaches them quickly and simply through every important step of the rescue process.

*\*Please call for more information.*

### Coolay Four Pegs – Bosworth

Made of special wood. Are shaped to fit into the occlusal grooves, permitting any restoration to be fully seated.

- |             |          |            |
|-------------|----------|------------|
| 50/Pkg..... | 135-0727 | 1@ \$26.40 |
|-------------|----------|------------|

### Brush-N-Floss – Hager Worldwide

*Study Model*

- Lightweight plastic, ivory and pink colours: weight only 320g.
- Twice the normal size with simple metal articulator
- For demonstration of tooth brushing and flossing, and patient motivation
- Model comes with extra large toothbrush.

- |            |          |             |
|------------|----------|-------------|
| Each ..... | 765-4767 | 1@ \$106.25 |
|------------|----------|-------------|



### Haricot Kidney Bassin – Almedic

Stainless Steel, 10"

- |            |          |            |
|------------|----------|------------|
| Each ..... | 105-2471 | 1@ \$21.95 |
|------------|----------|------------|



### HR® Binocular Loupes – HEINE

*With High Resolution technology for crisp, bright images with a super-wide field of view*

*\* Many configurations available.  
Please call for details*



### HRP® Binocular Loupes – HEINE

*With High Resolution Prismatic optics for high magnification and crystal clear images "edge-to-edge".*

*\* Many configurations available.  
Please call for details*



### F.O. LoupeLight® – HEINE

*The NEW Fiber Optics loupe Light from Heine.*

- Co-axial illumination and observation
- Light spot adjustable
- Weighs just 60g
- Water resistant housing

*\* Please call for more information*



### LED Lights – Pyramex

- Ideal for low light situations
- Fits hard hats and most safety glasses
- 2-bulb light with long battery life - over 40 hours

- |           |          |           |
|-----------|----------|-----------|
| Kit ..... | 246-3461 | 1@ \$8.00 |
|-----------|----------|-----------|

### Masque™ Silicone Lubricant – Bosworth

Non-soluble in the oral environment or monomer and not affected by head. Cannot run, tasteless, odourless, compatible with oral tissue, protects silicate fillings, lubricates rubber dams and forceps and can be applied whenever a lubricant or coating agent is required.

- |                  |          |            |
|------------------|----------|------------|
| 2.5oz. Tube..... | 138-4320 | 1@ \$23.45 |
|------------------|----------|------------|

### Mix Sticks – Zirc

*Disposable*

White plastic, double-end, with one square and one rounded tip.

- |              |          |            |
|--------------|----------|------------|
| 100/Pkg..... | 771-5972 | 1@ \$14.30 |
|--------------|----------|------------|

### Mini Lip Balm – Kamco

These mini lip balms are ideal for soothing dry chapped lips. Patients can use it while in the dental chair and then take it home for later use. They'll think of your office every time they apply it. Available in three options.

Give your patients a choice. Chap Ice Lip Balm comes in three great flavours: Orange, Kiwi Lime and Cherry



Chap Ice (assorted)  
Fishbowl of 100 .....218-6625      1@ \$75.50

An all time favourite, this all natural petroleum-free product contains refreshing Tea Tree Mint.



Lip Naturals  
Fishbowl of 100 .....218-6609      1@ \$75.50

A great way to reward your younger patients. Lip Rageous comes in four outrageous flavours: Strawberry, Watermelon, Red Apple and Blue Raspberry



Lip Rageous (assorted)  
Fishbowl of 100 .....218-6633      1@ \$75.50

### Seiler Loupes – Seiler

- Provides vision not attainable with the naked eye
- High resolution
- Optics are sealed to prevent dust and moisture
- Double hinge flip-up design
- All frames are 100% titanium
- Highly corrected optics that allow the entire tooth to be viewed
- Working distance: Short (340mm or 13.4"), Regular (420mm or 16.5") and Long (500mm or 19.7")

#### Each (Titanium)

Available in Black, Blue, Red and Silver

2.5X      1@ \$639.55  
3.0X      1@ \$704.20  
3.5X      1@ \$769.25



#### Each (Sport)

Available in Black, Blue and Silver

2.5X      1@ \$780.95  
3.0X      1@ \$846.55  
3.5X      1@ \$912.20



Please call for more details

www.pattersondental.ca

### Snoring Solution – MPowRx

- Invented by dentist & sleep researcher
- Clinically proven to be highly effective in reducing night-time apnea and snoring disturbances
- Health Canada cleared as Class I Medical Device
- No custom fitting required by dentist, making it easy to administer and support in your dental office
- Economical and accessible device to enable patients greater access to care
- Clinical study conducted by Dr. Leslie Dort (MSc, DDS, Dip ABDSM) and published in the Journal of Sleep and Breathing, May 2008"



Snoring Solution (6/Pkg) .....531-8993      1@ \$484.60  
Dental Starter Kit .....531-8985      1@ \$505.40  
Brochures (50/Pkg) ....531-8977      1@ \$19.25  
Brochure Holder .....531-8969      1@ \$6.25  
Laminated Poster .....531-8951      1@ \$19.25  
Posters & Stand.....531-8944      1@ \$12.50  
Device Cleaner (96/Pkg) .....531-8936      1@ \$36.95



### LED Loupe Light – Seiler

The Seiler LED Loupe Light is designed with precision engineering to provide the brightest light available.

#### Additional features are:

- Ultra Lightweight Design
- Constant light output as battery discharges
- Over 10,000 LUX
- Built-in battery gauge which displays the power level of the battery pack
- 7 intensity settings with up to 15 hours of illumination at the maximum setting
- Rechargeable Battery Pack
- Modular – will fit most flip-up hinge design Loupes

### Toothprints® – Kerr



#### Dental Identification Impression

An arch-shaped thermoplastic wafer that, when softened in hot water, registers a child's bite impression for identification and provides DNA and saliva for tracking them.

#### Introductory Kit

25/Box  
Child .....438-8104      1@ \$86.40



### TruRX Digital Denture Prescription – Dentsply

Dentsply TruRX, The Intelligent Denture System is an interactive program (patent pending) from Dentsply Prosthetics that standardizes and simplifies the denture consultation process, as it facilitates dental professional/patients and dental professional/laboratory communications. The premium program includes selections from our Trubyte® Bioform® moulds and IPN® material

TruRx Kit .....276-9008      1@ \$115.00  
Sheet TruRx Calibration  
Tools (set of 6) .....276-9016      1@ \$14.95



### Zekrya Gingival Protector – Zenith/DMG

#### Gingival Protector Kit

Zekrya Gingival Protector provides complete protection to your patient's gums during tooth preparation without sacrificing visibility or mobility. Two different size stainless steel arcs that shield the gingiva have been designed to fit the necks of all teeth. Contains: 1 swivel mandrel, 10 small arcs and 10 large arcs.

Kit .....324-2500      1@ \$296.45